



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
1	Achondroplasia	Diagnosis, FGFR3, Mutation Germline	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	Radiography/ medical impression/ reports	Y	Y
2	ALK Common Somatic Mutation	Somatic Mutation, Neuroblastoma	Mol- Trans	Sanger Sequencing	PBE or DNA or paraffin block	5-10 ml/> 20% Tumoral Tissue	2-8 °C	14	Pathology report/ Patient clinical information/ Form TE-REX-FO-097	N	N
3	Alpha Thal Female	Alpha Globin, Carrier Detection, Hematologic diseases, Alpha Thalassemia	Mol- Thalasemia	Sanger Sequencing- Multiplex PCR	PBE or DNA	5-10 ml	2-8 °C	30	CBC- Hb Electrophoresis-Ferritin	Y	N
4	Alpha Thal Male	Alpha Globin, Carrier Detection, Hematologic diseases, Alpha Thalassemia	Mol- Thalasemia	Sanger Sequencing- Multiplex PCR	PBE or DNA	5-10 ml	2-8 °C	30	CBC- Hb Electrophoresis-Ferritin	Y	N
5	Alpha Thal Other	Alpha Globin, Carrier Detection, Hematologic diseases, Alpha Thalassemia	Mol- Thalasemia	Sanger Sequencing- Multiplex PCR	PBE or DNA	5-10 ml	2-8 °C	30	CBC- Hb Electrophoresis-Ferritin	Y	N
6	Alpha-Beta-XMN	Alpha Globin, Beta- Globin, Thalassemia, Hematologic diseases, Thalassemia Diagnosis, Thalassemia Major	Mol- Thalasemia	Sanger Sequencing- ARMS PCR- Multiplex PCR	PBE or DNA	5-10 ml	2-8 °C	30	CBC- Hb Electrophoresis-Ferritin- Medical impression- Transfusion dependency	Y	N



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
7	ASXL1 Common Somatic Mutation	ASXL1, Somatic Mutation, Cancer, Myeloproliferative Neoplasms, (MPN)	Mol- Trans	Sanger Sequencing	PBE or DNA or paraffin block	5-10 ml/> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
8	ATRX Common Somatic Mutation	ATRX, Somatic Mutation, Nervous System Cancer	Mol- Trans	Sanger Sequencing	PBE or DNA or paraffin block	5-10 ml/> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
9	Beta Thal Female	Beta-Globin, Thalassemia; Hematologic diseases	Mol- Thalassemia	Sanger Sequencing- ARMS PCR	PBE or DNA	5-10 ml	2-8 °C	30	CBC- Hb Electrophoresis-Ferritin	Y	N
10	Beta Thal Male	Beta-Globin, Thalassemia; Hematologic diseases	Mol- Thalassemia	Sanger Sequencing- ARMS PCR	PBE or DNA	5-10 ml	2-8 °C	30	CBC- Hb Electrophoresis-Ferritin	Y	N
11	Beta Thal Other	Beta-Globin, Thalassemia; Hematologic diseases	Mol- Thalassemia	Sanger Sequencing- ARMS PCR	PBE or DNA	5-10 ml	2-8 °C	30	CBC- Hb Electrophoresis-Ferritin	Y	N
12	BRAF Mutation (Tissue)	BRAF, V600, Exon15, Somatic Mutation Nervous System Cancer, Melanoma, Colorectal Cancer, Lymphoma, Hairy cell leukemia	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE	5-10 ml/> 20% Tumoral Tissue	2-8 °C	7	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
13	CAH-Strip Assay	Diagnosis, big genitalia, CAH, Congenital adrenal hyperplasia	Mol- Others	Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
14	CALR	CALR, Somatic Mutation, Myeloproliferative Neoplasms (MPN), Cancer	Mol- Trans	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	Prescription, medical impression, CBC / reports	N	by the attending physician
15	Hereditary Diffuse Gastric Cancer (HDGC); CDH1	Gastric cancer, Sanger sequencing, Hereditary cancer	Mol-Gene seq	Sanger Sequencing	PBE	10 ml	2-8 °C	60	medical impression/ reports	Y	Y
16	CEBPA Common Somatic Mutation	CEBPA, Somatic Mutation, Acute Myeloid Leukemia (AML), Cancer	Mol- Trans	Sanger Sequencing	PBE or DNA or paraffin block	5-10 ml/> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
17	CF Common Mutation Female	CFTR gene (Cystic Fibrosis)	Mol- Others	Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
19	CF Common Mutation Other	CFTR gene (Cystic Fibrosis)	Mol- Others	Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک موردنیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
20	CKIT	CKIT, Somatic Mutation, Acute Myeloid Leukemia (AML), Myeloproliferative Neoplasms (MPN),Sarcoma, Melanoma, Cancer	Mol- Trans	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
21	CKIT-Tissue	CKIT, Somatic Mutation, Gastric Cancer, GIST , Cancer	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
22	Clotting Panel	MTHFR (A1298C) (C677T), PAI, FV;FXIII(V34L in Factor XIII gene), Prothrombin (Factor II) (G20210A),	Mol- Others	Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	by the attending physician
23	Connexin 26 Female	GJB2- Connexin26-35Del G, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Audiometry result	Y	Y
24	Connexin 26 Male	GJB2- Connexin26-35Del G, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Audiometry result	Y	Y
25	Connexin 26 Other	GJB2- Connexin26-35Del G, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Audiometry result	Y	Y
26	Connexin 30	GJB6- Connexin30, Diagnosis, Carrier	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Audiometry result	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
27	CYP2C9 Genotype	MAYZENT, Multiple Sclerosis, MS	Mol- Trans	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	Y	N
28	DM1	DM1	Mol- Others	Fragment Analysis	PBE or DNA	5-10 ml	2-8 °C	30	medical impression/ reports	Y	Y
29	DNA Bank	DNA Bank	Mol- Extraction	DNA Extraction	PBE or Tis	5-10 ml Blood or 20% Tumoral Tissue	2-8 °C	1095	Prescription	Y	Y
30	DNA Collection- mouth swab		Mol- Extraction	DNA Extraction			2-8 °C			Y	N
31	DNA-Extraction- Report	Checking of Quality and Quantity of DNA	Mol- Extraction	DNA Extraction	PBE or Tis	5-10 ml Blood or 20% Tumoral Tissue	2-8 °C	3	medical impression/ reports	Y	N
32	DNMT3A Common Somatic Mutation	DNMT3A, Common Somatic Mutation, Pharmacogenetics, Myeloproliferative Neoplasms (MPN), Acute Myeloi Leukemia (AML)	Mol- Trans	Sanger Sequencing	PBE or Tis or DNA or FFPE	5-10 ml/> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
33	DPYD (5-fu)	colorectal cancer , Solid tumor , IVS14+1G>A	Mol- Others	Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	by the attending physician



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
34	Duchenne-Ind	Duchenne Muscular Dystrophy, DMD, BMD, Dystrophin. Deletion /Duplication, Diagnosis, Carrier Detection	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
35	EGFR	Lung Cancer, Non Small Cell Carcinoma, Somatic Mutation, Exon18, 19, 20, 21	Mol- Trans	Sanger Sequencing/ Reverse Strip Assay	Tis or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	7	medical impression/ reports/Pathology Reports/prescription	N	N
36	ESR1	Breast cancer, Prostat cancer, Estrogen receptor	Mol- Trans	Sanger Sequencing	Tis or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	N
37	Factor XIII	FXIII(V34L in Factor XIII gene, Diagnosis, Carrier Detection	Mol- Others	Sanger Sequencing/ Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	N
38	Factor-V (Leiden)	Factor-V (Leiden) - G1691A	Mol- Others	Sanger Sequencing/ Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	N
39	FancA-Seq-MLPA	FancA. Fanconi anemi, Germline Mutation	Mol- Gene Sequencing	Sanger Sequencing/MLPA	PBE or DNA	5-10 ml	2-8 °C	60	medical impression/ reports	Y	Y
40	FAP (APC gene mutation analysis)	Cancer, Pericardial, Peritoneal, Screening, Slide, Urine	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	60	-	Y	by the attending physician



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
41	FLT-3	FLT3-Acute Myeloid Leukemia (AML), FLT3-ITD, FLT3-835-Somatic Mutation	Mol- Trans	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	CBC/medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
42	FMF	FMF-Familial Mediterranean Fever, MeFv Mutation, Common Mutation, Diagnosis, Carrier Detection	Mol- Others	Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	Y
43	Fragment Reaction	Diagnostic test/ Research	Mol- Sequencing	Fragment Analysis	PCR Product	25UL	2-8 °C	3	لیبل نمونه روی میکروتیوپ و اطلاعات فرگمنت در اکسل ارسال	N	N
44	Fragment Run	Diagnostic test/ Research	Mol- Sequencing	Fragment Analysis	PCR Product	25UL	2-8 °C	3	لیبل نمونه روی میکروتیوپ و اطلاعات فرگمنت در اکسل ارسال	N	N
45	FraX-Mol	FraX, Premature Ovarian Failure (POF), X link Syndromic Mental Retardation, Diadnosis, Carier Detection	Mol- Others	PCR- Fragment Analysis/Southern	PBE or DNA	5-10 ml	2-8 °C	30	medical impression/ reports	Y	Y
46	Friedreich`s ataxia	FXN, GAA Expanation, Friedreich`s ataxia	Mol- Others	PCR- Fragment Analysis	PBE or DNA	5-10 ml	2-8 °C	30	medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
47	GATA Common Somatic Mutation	Common Somatic Variant, Myeloproliferative	Mol- Trans	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	Y
48	Gaucher	GBA gene;beta-glucocerebrosidase, Dagnosis, Carrier Detection	Mol- Others	Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	Y
49	Gene Seq 02-05 Exon	Gemline Mutation, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Prescription, medical impression/ reports	Y	Y
50	Gene Seq 06-10 Exon	Gemline Mutation, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	medical impression/ reports	Y	Y
51	Gene Seq 11-25 Exon	Gemline Mutation, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	60	medical impression/ reports	Y	Y
52	Gene Seq 26-50 Exon	Gemline Mutation, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	60	medical impression/ reports	Y	Y
53	H3F3A Common Somatic Mutation	H3F3A, Common Somatic Mutation, Pharmacogenetics, Nervous System Cancer, Cancer	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
54	Hemochromatosis	HFE, Diagnosis, Carrier Detection	Mol- Others	Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
55	Hemophilia A- Inversion 22	Inversion 22	Mol- Gene Sequencing	Long Range PCR	PBE or DNA	5-10 ml	2-8 °C	30	medical impression/ reports/Pathology Reports/prescription	Y	Y
56	Hemophilia A- Sequencing	Sequencing	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	60	medical impression/ reports/Pathology Reports/prescription	Y	Y
57	Huntington	HTT gene, CAG Expanation, Huntington's disease (HD), autosomal dominant inheritanc	Mol- Others	PCR- Fragment Analysis	PBE or DNA	5-10 ml	2-8 °C	30	medical impression/ reports	Y	Y
58	Hypochondroplassi a	FGFR3, Germline Mutation, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
59	IDH1 Common Somatic Mutation	IDH1, Somatic Mutation, Nervous system Cancer, Myeloprolifrative Neoplasms (MPN)	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
60	IDH2 Common Somatic Mutation	IDH2, Somatic Mutation, Nervous system Cancer, Myeloprolifrative Neoplasms (MPN)	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
61	Imatinib Resistance	CML, T315I, resistance to tyrosine kinase inhibitor (TKI) therapy,	Mol- Trans	Sanger Sequencing	PBE or cDNA	5-10 ml	2-8 °C	14	CBC/ Prescription/ medical impression/ reports	N	by the attending physician



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
62	Irinotecan Toxicity	UGT1A1*28 genotype, Pharmacogenomics, Pancreatic Neoplasms, Colorectal Neoplasms,	Mol- Others	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	جواب غیر طبیعی
63	JAK2 Exon 12	JAK2, Myeloproliferative Neoplasms (MPN), Leukemia, MPN	Mol- Trans	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	CBC/ Prescription/ medical impression/ reports	N	by the attending physician
64	JAK2 gene v617F Mutation	JAK2, V617F, Leukemia, Myeloproliferative Neoplasms (MPN),CML,AML	Mol- Trans	Sanger Sequencing/ Real time PCR	PBE or DNA	5-10 ml	2-8 °C	7	CBC/ Prescription/ medical impression/ reports	N	by the attending physician
65	K-RAS	KRAS Common Mutation, Exon 2, 3, 4, Codone 12-13.61- colorectal Cancer	Mol- Trans	Sanger Sequencing/ Reverse Strip Assay/Real time PCR	Tis Or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	7-10	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
66	MCC+QF-PCR Aneuploidies	Maternal Cell Contamination	Mol- Trans	Fragment Analysis	AF/AFC/Tis/P BE or DNA	5-10 ml	2-8 °C	3	Prescription/ medical impression/ reports	Y	Y
67	MCC-PND	Maternal Cell Contamination	Mol- Trans	Fragment Analysis	PBE or DNA	5-10 ml	2-8 °C	3	Prescription/ medical impression/ reports	Y	Y
68	MLPA- CAPN3	MLPA- Large Deletion /Insertion, Diagnosis, Carrier Detection, CAPN3, CALpain3	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
69	MLPA- Celiac DQ2 & DQ8	Celiac disease (HLA- DQ2/ HLA-DQ8)	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	Prescription/ medical impression/ reports	Y	N
70	MLPA-Alpha	MLPA, Alpha Thalassemia, Detection of copy number changes in the Alpha globin gene cluster	Mol- Thalasemia	MLPA	PBE or DNA or CV or AF or AFC	5-10 ml	2-8 °C	14	CBC- Hb Electrophoresis	Y	Y
71	MLPA- Aneuploidies	MLPA-Aneuploidies, QF, Down Syndrome	Mol- Others	MLPA	PBE or DNA or CV or AF or AFC	5-10 ml	2-8 °C	7	Prescription/ medical impression/ reports	Y	Y
72	MLPA-Beta	MLPA, Delta-Beta Thalassemia, Detection of copy number changes in the Beta globin gene cluster	Mol- Thalasemia	MLPA	PBE or DNA or CV or AF or AFC	5-10 ml	2-8 °C	14	CBC- Hb Electrophoresis	Y	Y
73	MLPA-BRCA I-II	MLPA, Breast and Ovarian Cancer, Diagnosis, Carrier Detection	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
74	MLPA-CLL	deletion, Insertion, Chronic lymphocytic leukemia	Mol- Trans	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	Prescription/ medical impression/ reports	Y	Y
75	MLPA-CMT (Duplication- PMP22)	MLPA, Charcot-Marie-Tooth disease (CMT), HNPP, Diagnosis, Carrier Detection	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک موردنیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
76	MLPA-Del/Dup I	MLPA- Large Deletion /Insertion, Diagnosis, Carrier Detection	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	120	medical impression/ reports	Y	Y
77	MLPA-MGMT Methylation	Hypermethylation, glioma	Mol- Trans	MLPA	Tis Or DNA or FFPE	> 20% Tumoral Tissue/5-10 ml	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
78	MLPA- Microdeletions	MLPA- Deletion /Insertion, Diagnosis, Carrier Detection 1p36 deletion syndrome/2p16.1-p15 microdeletion syndrome/2q23.1 microdeletion/microduplic ation syndrome/Glass syndrome/3q29 microdeletion syndrome/3q29 microduplication syndrome/Wolf- Hirschhorn syndrome/Cri- du-Chat syndrome/Sotos syndrome/Williams- Beuren syndrome/Williams- Beuren duplication syndrome/Langer- Giedion syndrome/9q22.3 microdeletion syndrome/DiGeorge syndrome-2/Prader-Willi syndrome/Angelman syndrome/Witteveen- Kolk* / 15q24 microdeletion syndrome/Rubinstein- Taybi syndrome/Miller	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	N



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
79	MLPA-MLH1-2,PMS2,MSH6 Hypermethylation	Hypermethylation, Somatic Mutation, Mis match repair gene, microsatellitc instability	Mol- Trans	MLPA	Tis Or DNA or FFPE	> 20% Tumoral Tissue/5-10 ml	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	Y	by the attending physician
80	MLPA-Mutation CAH	Diagnosis, big genitalia, CAH, Congenital adrenal hyperplasia	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
81	MLPA-SGCA	Limb-girdle muscular dystrophy; MLPA- Large Deletion Duplication, Diagnosis, Carrier Detection	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
82	MPL-Mutation	MPL, Somatic Mutation- Myeloproliforative Neoplasms (MPN), Cancer	Mol- Trans	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	CBC/ Prescription/ medical impression/ reports	N	by the attending physician
83	MSI	Micro Satellite Inastability- Colorectal Cancer , MMR, Miss Match Repare, HNPCC	Mol- Trans	Fragment Analysis	PBE and Tis or FFPE	5-10 ml/> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
84	MTH-FR	MTHFR (A1298C) (C677T)	Mol- Others	Sanger Sequencing/ Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	by the attending physician



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک موردنیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
85	Mutation Confirmation-Family (up to 3 person)	Mutation Confirmation, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
86	Mutation Confirmation-Ind (additional person)	Mutation Confirmation, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
87	NF1 Common Somatic Mutation	NF1, Somatic Mutation, Sarcoma, Cancer	Mol- Trans	Sanger Sequencing	Tis Or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
88	NGS- Carrier Screening-Medical Exome	Carrier Screening, OMIM Genes	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports/pedigree	Y	Y
89	NGS- Multi Cancer With CNV Panel (Solid Tumor and Non solid Tumor Cancer)	Solid Tumor, Cancer, NGS, Germline, Somatic	Mol- NGS	NGS	PBE or DNA or Tis or FFPE	5-10 ml/> 20% Tumoral Tissue	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
90	NGS-Autism Spectrum Disorders-Exome	Autism, ASD, intellectual disability	Mol- NGS	NGS	PBE or DNA or Tis or FFPE	5-10 ml/> 20% Tumoral Tissue	2-8 °C	120	Prescription/ medical impression/ reports	Y	Y
91	NGS-Ataxia Exome	NGS, Ataxia, Germline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
92	NGS-BRCA I-II Panel	NGS- Germline mutation, Diagnosis, Breast cancer, Ovarian Cancer, pancreatic cancer, BRCA1- BRCAII, BRCA1, BRCA2, Sequencing	Mol- Gene Sequencing	NGS	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
93	NGS-BRCA I-II Plus del/dup	NGS- Germline mutation, Diagnosis, Breast cancer, Ovarian Cancer, pancreatic cancer, BRCA1- BRCAII, BRCA1, BRCA2, Sequencing	Mol- Gene Sequencing	NGS	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
94	NGS-BRCA I-II Plus del/dup (Somatic) Panel	NGS- Germline mutation, Diagnosis, Breast cancer, Ovarian Cancer, pancreatic cancer, BRCA1- BRCAII, BRCA1, BRCA2, Sequencing	Mol- Trans	NGS	DNA in streck tubes	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
95	NGS-Carrier Screening Female	NGS,Parents (Including carrier detection), Gemline Mutation, Diagnosis, Carrier Detection, Reproductive Screening	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
96	NGS-Carrier Screening Female (Without DMD)	NGS,Parents (Including carrier detection), Gemline Mutation, Diagnosis, Carrier Detection, Reproductive Screening	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
97	NGS-Carrier Screening male	NGS,Parents (Including carrier detection), Gemline Mutation, Diagnosis, Carrier Detection, Reproductive Screening	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
98	NGS-Ciliopathies Exome	NGS-Ciliopathies, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
99	NGS-CMT & Sensory Neuropathies (103G) Exome	NGS, Sensory Neuropathies, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
100	NGS-CNV Exome	NGS, Gemline Mutation, Diagnosis, Carrier Detection, Solid Tumor, Cancer	Mol- NGS	NGS	PBE or DNA or Tis or FFPE	5-10 ml/> 20% Tumoral Tissue	2-8 °C	120	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
101	NGS-Corpus Callosum Agenesis Exome	NGS-Corpus Callosum Agenesis, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
102	NGS-Craniofacial anomalies Exome	NGS, Carniofacial Anomalies, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
103	NGS-Disorders in sex development Exome	NGS, Disorders in sex development, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
104	NGS-Duo (Proband & 1 family member) Exome	NGS, Duo (Proband & 1 family member, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
105	NGS-Epilepcy Exome	NGS-Epilepcy, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
106	NGS-Exome (Medical exome for OMIM refrenced genes 200X)	NGS,Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
107	NGS-Exome (Medical exome for OMIM referenced genes)	NGS-Exome (Medical exome for OMIM referenced genes), Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
108	NGS-Exome+ Mito 100X	Mitochondrial disorder, heteroplasmy, mt DNA Screening	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
109	NGS-Hearing Loss Exome	NGS-Hearing Loss, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
110	NGS- Hemostatic/Thro mbotic disorders Exome	NGS- Hemostatic/Thrombotic disorders, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
111	NGS-Hereditary Cancer	NGS, Hereditary Cancer, Gemline Mutation, Diagnosis, Carrier Detection	Mol-Gene seq	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
112	NGS-Hereditary comprehensive cancer panel	NGS, Hereditary Cancer,	Mol-Gene seq	NGS	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک موردنیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
113	NGS-Hereditary Cancer plus CNV	NGS,Hereditary Cancer, Gemline Mutation, Diagnosis, Carrier Detection,CNV	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
114	NGS-Hereditary Spastic Paraplegia Exome	NGS, Hereditary Spastic Paraplegia Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
115	NGS-ID (Syndromic-NonSyndromic) Exome	NGS-ID (Syndromic-NonSyndromic), Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
116	NGS-IBMFS (Aplastic anemia)	Aplastic anemia, Fanconi anemia, shmachman diamond syndrome, diamond blackfan syndrome, bone marrow failure	Mol-Gene seq	NGS	PBE or DNA or Tis or FFPE	5-10 ml/> 20% Tumoral Tissue	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
117	NGS-Iron disorders Exome	NGS,Iron disorders, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک موردنیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
118	NGS-Mendelian Disorders Exome	NGS, Mendelian Disorders, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
119	NGS-Metabolic disorders Exome	NGS-Metabolic disorders, Metabolic disorders	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
120	NGS- Mitochondrial disorders Exome	NGS,Mitochondrial disorders, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
121	NGS-Movement disorders Exome	NGS,Movement disorders, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
122	NGS-Multiple congenital anomalies Exome	NGS,Multiple congenital anomalies, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
123	NGS-Muscular Dystrophy & Muscledisorders Exome	NGS,Muscular Dystrophy & Muscledisorders, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
124	NGS-Onco Myeloid Panel	KDM6A, Leukemia, NGS,Cancer,TP53,SMC1 A,PTEN,NMP1,HRAS,ETV6,CEBPA,BRAF,FLT3,JAK2,U2AF1,WT1,SF3B1,ATRX,ASXL1,K1T,NRAS,E2H2,IDH1,IDH2	Mol-Trans	NGS	PBE or DNA	10 ml	2-8 °C	30	Prescription/ medical impression/ reports	Y	Y
125	NGS-Onco Multi cancer RNA fusion panel	Leukemia,RNA,NGS,Cancer,RAF1,BRAF,NTRK1,2,3,EGFR,FGFR3,PBX1,NRG,ROS1,RET,ALK,TFE3,PRKACA,ERG,PPARG,MET	Mol-Trans	NGS	PBE or FFPE	10 ml	2-8 °C	30	Prescription/ medical impression/Pathology reports	Y	Y
126	NGS-Parents (Including carrier detection) Exome	NGS,Parents (Including carrier detection), Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
127	NGS-Primary immunodeficiencies Exome	NGS-Primary immunodeficiencies, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
128	NGS-Reanalysis of data	NGS,Reanalysis of data, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
129	NGS-Research Exome 100X	NGS,Research, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
130	NGS-RP Exome	NGS,Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
131	NGS-Skin disorders Exome	NGS,Skin Disorders, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
132	NGS-TMB and MSI	Thyroid cancer, Lung, Colorectal cancer, EGFR Mutation	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	45	Prescription/ medical impression/ reports	Y	Y
133	NGS-Trio (Proband & 2 family preferred members) Exome	NGS-Trio (Proband & 2 family preferred members) Exome, Gemline Mutation, Diagnosis, Carrier Detection	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
134	NGS-Vision disorders Exome	NGS,Gemline Mutation, Diagnosis, Carrier Detection, Vision disorders	Mol- NGS	NGS	PBE or DNA	5-10 ml	2-8 °C	90-120	Prescription/ medical impression/ reports	Y	Y
135	NGS-Whole Genome Sequencing (WGS)	Whole Genome Sequencing, Breakpoint detection	Mol- NGS	Whole genome sequencing	PBE or DNA	5-10 ml	2-8 °C	90	medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
136	NPM-1	NPM1, Somatic Mutation	Mol- Trans	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	7	Prescription/medical impression/ reports/CBC	N	by the attending physician
137	N-RAS	NRAS Common Mutation, Exon 2, 3, 4, Codone 12-13.61- colorectal Cancer	Mol- Trans	Sanger Sequencing/ Reverse Strip Assay/Real time PCR	Tis Or DNA	> 20% Tumoral Tissue	2-8 °C	7-10	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
138	PAI	PAI, Clotting Problem, Trombosis	Mol- Others	Sanger Sequencing/ Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	by the attending physician
139	Panel ALL	t(1;19) - t(12, 21) - t(9;22); t(4;11), Translocation , Leukemia, Diagnosis, Prognosis	Mol- Trans	RT-PCR/Real Time PCR	PBE or BM	5-10 ml	2-8 °C	14	All disease diagnosis documents - prescription - CBC	N	by the attending physician
140	Panel AML	t(4;11) - t(8, 21) - t(9;22); t(15;71), inv16, Translocation , Leukemia, Diagnosis, Prognosis	Mol- Trans	RT-PCR/Real Time PCR	PBE or BM	5-10 ml	2-8 °C	14	All disease diagnosis documents - prescription - CBC	N	by the attending physician
141	Paternity	Paternity	Mol- Trans	Fragment Analysis	PBE or DNA	5-10 ml	2-8 °C	14	Identity documents	Y	Y
142	Paternity-Ind	Paternity	Mol- Trans	Fragment Analysis	PBE or DNA	5-10 ml	2-8 °C	14	Identity documents	Y	Y
143	PCR- Ampilification	PCR	Mol- Trans	PCR	DNA in streck tubes	5-10 ml	2-8 °C	3	Prescription/ medical impression/ reports	N	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
144	PDGFRa-Tissue	PDGFRa, Somatic Mutation, GIST (Exon18), Sarcoma	Mol- Trans	Sanger Sequencing	Tis Or DNA	> 20% Tumoral Tissue	2-8 °C	14	medical impression/ reports/Pathology Reports/prescription	N	by the attending physician
145	PGD	PGD-PGT-M-Monogenic disease-IVF-PND	Mol- PGD	Sanger sequencing/ Linkage analysis	Blastomere/ blastocyst	1-4 Blastomeres	2-8 °C	2	medical impression/ reports	Y	Y
146	PGD I	PGD-PGT-M-Monogenic disease-IVF-PND	Mol- PGD	Sanger sequencing/ Linkage analysis	Blastomere/ blastocyst	<4 Blastomeres	2-8 °C	2	medical impression/ reports	Y	Y
147	PGD II	PGD-PGT-M-Monogenic disease-IVF-PND	Mol- PGD	Sanger sequencing/ Linkage analysis	Blastomere/ blastocyst	>4 Blastomeres	2-8 °C	2	medical impression/ reports	Y	Y
148	PGD SEX I	PGD-PGT-M-Monogenic disease-IVF-PND-Pre-PGT-M-SEX determination	Mol- PGD	AMX/SRY amplification-QF analysis	Blastomere/ blastocyst	1-4 Blastomeres	2-8 °C	2	medical impression/ reports	Y	Y
149	PGD SEX II	PGD-PGT-M-Monogenic disease-IVF-PND-Pre-PGT-M-SEX determination	Mol- PGD	AMX/SRY amplification-QF analysis	Blastomere/ blastocyst	>4 Blastomeres	2-8 °C	2	medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک موردنیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
150	PGD-Pre	PGD-PGT-M-Monogenic disease-IVF-PND-Pre- PGT-M	Mol- PGD	Sanger sequencing/ Linkage analysis	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
151	PGD-Pre Beta & GJB2	PGD-PGT-M-Monogenic disease-IVF-PND-Pre- PGT-M	Mol- PGD	Sanger sequencing/ Linkage analysis	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
152	PKU Female	PKU (Phenylketonuria), PAH Gene, Germline Mutation, Diagnosis, Carrier Detection, Metabolic disorders	Mol- Gene Sequencing	Linkage analysis	PBE or DNA	5-10 ml	2-8 °C	30	Prescription/ medical impression/ reports	Y	Y
153	PKU Male	PKU (Phenylketonuria), PAH Gene, Germline Mutation, Diagnosis, Carrier Detection, Metabolic disorders	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Prescription/ medical impression/ reports	Y	Y
154	PKU Other	PKU (Phenylketonuria), PAH Gene, Germline Mutation, Diagnosis, Carrier Detection, Metabolic disorders	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
155	PND- PGD	PGD-PGT-M-Monogenic disease-IVF-PND-Pre-PGT-MPreimplantation genetic diagnosis	Mol- PGD	Sanger Sequencing	CVS/Amnion	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
156	PND-Alpha	PND- Alpha Thalassemia, Prenatal Diagnosis	Mol- Thalasemia	Sanger Sequencing- Multiplex PCR- Alpha MLPA	PBE or DNA or CV or AF or AFC	5-10 ml	2-8 °C	14	CBC- Hb Electrophoresis-Ferritin	Y	Y
157	PND-CF	PND -CF (Cystic fibrosis), Sequencing, MLPA,Prenatal Diagnosis	Mol- Others	Sanger Sequencing/ Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
158	PND-Duchenne	PND -DMD , Sequencing, MLPA, Prenatal Diagnosis	Mol- Others	Sanger Sequencing/ML PA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
159	PND-FraX	PND -FRAx , PreMutation, Prenatal Diagnosis	Mol- Others	Fragment Analysis/Southern	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
160	PND-Gaucher	PND -Gaucher Germline Mutation, Prenatal Diagnosis	Mol- Others	Sanger Sequencing/ Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y
161	PND-Mutation Confirmation	Germline Mutation, Sequencing, Prenatal Diagnosis	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	Prescription/ medical impression/ reports	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
162	PND-Others	PND, Prenatal Diagnosis	Mol- Others	Sanger Sequencing-MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/reports	Y	Y
163	PND-PKU	PND-PKU, (Phenylketonuria), Germline Mutation, Prenatal Diagnosis, Metabolic disorders	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	Prescription/ medical impression/ reports	Y	Y
164	PND-PRE	PND, NEW Variant	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	60	Complete information of previous tests/clinical records/paraclinic tests	Y	Y
165	PND-Sex	PND- Sex Determination, Prenatal Diagnosis	Mol- Others	PCR	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/reports	Y	Y
166	PND-Sickle	PND- Sickle, Prenatal Diagnosis	Mol- Thalasemia	Sanger Sequencing	PBE or DNA or CV or AF or AFC	5-10 ml	2-8 °C	14	CBC- Electrophoresis- Ferritin	Y	Y
167	PND-SMA	PND, SMA, Prenatal Diagnosis	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/reports	Y	Y
168	PND-Thal Beta	PND, Beta thalassemia, Prenatal Diagnosis	Mol- Thalasemia	Sanger Sequencing- Beta MLPA	PBE or DNA or CV or AF or AFC	5-10 ml	2-8 °C	14	CBC- Electrophoresis- Ferritin	Y	Y
169	POLE Common Somatic Mutation	Somatic Mutatin, Endometrial Cancer POLE, Helisomatic mutations	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical compression/reports/pathology results/form TE-REX-FO-097/	Y	Y



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک موردنیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
170	PT-Pro- Thrombin (Factor II)	Factor II,Prothrombin,PT	Mol- Others	Sanger Sequencing/ Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	by the attending physician
171	QF-PCR Aneuploidies	Determination of Aneuploidy of Chr. 13, 18, 21 , x and Y in fetus, abnormality, Screening, Aneuploidies	Mol- Trans	Fragment Analysis	PBE or DNA and AF or PBE or DNA and CVs or PBE or DNA and TC	3-5 ml Amnion or CVS and 5- 10 ml PBE	2-8 °C	3	Prescription/ medical impression/ reports	Y	Y
172	RET Common Somatic mutation	Common Somatic Variant, Thyroid medullary Carcinoma, proto- oncogen, receptor tyrosin kinase	Mol- Trans	Sanger Sequencing	PBE or Tis or DNA	5-10 ml/> 20% Tumoral Tissue	2-8 °C	14	medical compression/reports/path ology results/prescription	N	by the attending physician
173	RT-PCR (1-19)	RT-PCR- ALL , Leukemia, Diagnosis	Mol- Trans	RT-PCR	PBE or BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
174	RT-PCR (12-21)	RT-PCR- ALL , Leukemia, Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE or BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
175	RT-PCR (15-17)	RT-PCR- AML, AML- M3, PML-RARA , Leukemia, Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE or BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
176	RT-PCR (15-17) Quantitative PMLRARA	RT-PCR- AML, AML- M3, PML-RARA , Leukemia, Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE or BM	5-10 ml	2-8 °C	14	All disease diagnosis documents/ prescription / CBC	N	by the attending physician



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
177	RT-PCR (4-11)	RT-PCR- ALL , Leukemia, Diagnosis	Mol- Trans	RT-PCR	PBE or BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
178	RT-PCR (8-21)	RT-PCR- AML , Leukemia,Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE or BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
179	RT-PCR (9-22) (P190KD-e1a2)	RT-PCR- CML, P190 BCR-ABL , Leukemia, Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE- BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
180	RT-PCR (9-22) (P190KD- P210KD- P230KD)	RT-PCR- CML, P190 BCR-ABL , Leukemia, Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE- BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
181	RT-PCR (9-22) (P210KD-b2a2- b3a2)	RT-PCR- CML, P210 BCR-ABL , Leukemia, Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE- BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
182	RT-PCR (9-22) (P230KD-e19a2)	RT-PCR- CML, P210 BCR-ABL , Leukemia, Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE- BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
183	RT-PCR (9-22) Quantitative (P190KD-e1a2)	RT-PCR- CML, P190 BCR-ABL , Leukemia, Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE- BM	5-10 ml	2-8 °C	14	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
184	RT-PCR (9-22) Quantitative (P210KD-b2a2- b3a2)	RT-PCR- CML, P210 BCR-ABL , Leukemia, Diagnosis	Mol- Trans	RT-PCR/Real Time PCR	PBE- BM	5-10 ml	2-8 °C	14	All disease diagnosis documents/ prescription / CBC	N	by the attending physician



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
185	RT-PCR inv (16)	RT-PCR- AML , Leukemia, Diagnosis	Mol- Trans	RT-PCR	PBE- BM	5-10 ml	2-8 °C	7	All disease diagnosis documents/ prescription / CBC	N	by the attending physician
186	Sequencing CAPN3	Sequencing, limb-girdle muscular dystrophy-1	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Prescription/ medical impression/ reports	Y	Y
187	Sequencing CF	CFTR-Cyctic Fibrosis- Germline Mutation, Sequencing, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
188	Sequencing DMD/BMD	DMD, Duchenne muscular dystrophy - Germline Mutation, Sequencing, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	60	Prescription/ medical impression/ reports	Y	Y
189	Sequencing FMF	MEFV, FMF- Germline Mutation, Sequencing, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Prescription/ medical impression/ reports	Y	Y
190	Sequencing Hematochromatosis	Germline Mutation, Sequencing, Diagnosis, Carrier Detection	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Prescription/ medical impression/ reports	Y	Y
191	Sequencing Reaction	Sequencing RUN	Mol- Sequencing	Sanger Sequencing	PCR Product	25UL	2-8 °C	3	لیبل نمونه روی میکروتیوپ، ارسال پرایمر و اطلاعات ژن و پرایمر در اکسل ارسالی	N	N



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی(روز) TAT (Day)	مدارک موردنیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
192	Sequencing Run	Sequencing Run	Mol- Sequencing	Sanger Sequencing	PCR Product	25UL	2-8 °C	3	لیبل نمونه روی میکروتیوپ، ارسال پرایمر و اطلاعات ژن و پرایمر در اکسل ارسالی	N	N
193	Sequencing TP53	Li-fraumeni Syndrome, P53, TP53, Germline Mutation, Sequencing, Diagnosis, Tumer protein P53, TRP53	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	30	Prescription/ medical impression/ reports	Y	Y
194	SEX- Determination	Sex-Determination	Mol- Others	PCR	PBE or DNA	5-10 ml	2-8 °C	5	medical impression/ reports	Y	Y
195	SF3B1 Common Somatic Mutation	Cancer,Sanger Sequencing, MDS, جهش های هلی سوماتیک، Breast Cancer, CLL, Endometrial Cancer	Mol- Trans	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	N	Y
196	Sickle Female	Sikle cell Thalassemia	Mol- Thalasemia	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	CBC- Hb Electrophoresis-Ferritin	Y	N
197	Sickle Male	Sikle cell Thalassemia	Mol- Thalasemia	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	CBC- Hb Electrophoresis-Ferritin	Y	N
198	Sickle Other	Sikle cell Thalassemia	Mol- Thalasemia	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	CBC- Hb Electrophoresis-Ferritin	Y	N



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
199	SLOS Gene Sequencing	Smith Lemli Opitz Syndrome, deficiency in the enzyme 7-dehydrocholesterol reductase, abnormality in cholesterol metabolism	Mol- Gene Sequencing	Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	Prescription/ medical impression/ reports	Y	جواب غیر طبیعی
200	SMA Female	SMA-Spinal Muscular Atrophy, Diagnosis, Carrier Detection	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	جواب غیر طبیعی
201	SMA FraX-Carrier screening female	Reproductive Screening, Pre-Pregnancy, Germline Screening, Carrier Screening	Mol- Others	MLPA, TP-PCR	PBE or DNA	5-10 ml	2-8 °C	30	Prescription/ medical impression/ reports	Y	جواب غیر طبیعی
202	SMA Male	SMA-Spinal Muscular Atrophy, Diagnosis, Carrier Detection	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	جواب غیر طبیعی
203	SMA Other	SMA-Spinal Muscular Atrophy, Diagnosis, Carrier Detection	Mol- Others	MLPA	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	جواب غیر طبیعی
204	SRSF2 Common Somatic Mutation	SRSF2, Somatic Mutation, Myeloproliferative Neoplasms (MPN), Cancer, Pharmacogenetics	Mol- Trans	Sanger Sequencing	PBE or Tis or DNA or FFPE	5-10 ml/> 20% Tumoral Tissue	2-8 °C	14	medical compression/ reports/ pathology results/ doctor's prescription	Y	by the attending physician
205	SRY/Sex Determination	SRY, Sex Determination	Mol- Others	PCR	PBE or DNA	5-10 ml	2-8 °C	3	medical impression/ reports	Y	جواب غیر طبیعی



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
206	TERT Common Somatic Mutation	TERT, Somatic Mutation, Nervous System Cancer, Cancer, Pharmacogenetics	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical compression/ reports/ pathology results/ doctor's prescription	Y	by the attending physician
207	Thal II-1-XMN	Thalassemia, XMN, Carrier Detection, Thalassemia Major	Mol- Thalasemia	ARMS PCR- XMNI Sanger Sequencing	PBE or DNA	5-10 ml	2-8 °C	14	CBC- Hb Electrophoresis	Y	N
208	TP53 Common Somatic Mutation	TP53, Somatic Mutation, Chronic Lymphocytic Leukemia (CLL), Myeloproliferative Neoplasms, Cancer, Pharmacogenetics	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical compression/ reports/ pathology results/ doctor's prescription	N	by the attending physician
209	TPMT	thiopurine therapy, thiopurine methyltransferase	Mol- Trans	Sanger Sequencing	PBE	5-10 ml	2-8 °C	14	medical impression/ reports	N	by the attending physician
210	U2AF1 Common Somatic Mutation	U2AF1, Common Somatic Mutation, Pharmacogenetics,, Myeloproliferative Neoplasms (MPN), Cancer	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE	> 20% Tumoral Tissue	2-8 °C	14	medical compression/ reports/ pathology results/ doctor's prescription	N	by the attending physician
211	Warfarin Resistance	Warfarin Resistance	Mol- Others	Reverse Strip Assay	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	by the attending physician



ردیف	Test	کلیدواژه keyword	ID	روش انجام method	نوع نمونه sample	حجم نمونه volume	شرایط نگهداری storage	زمان جوابدهی (روز) TAT (Day)	مدارک مورد نیاز document	نیاز به مشاوره قبل از انجام آزمایش دارد؟ Pre test counseling	نیاز به مشاوره بعد از انجام آزمایش دارد؟ Post test counseling
212	WGA PGD I		Mol- PGD	amplification	Blastomere/ blastocyst	1-4 Blastomeres	2-8 °C	2	medical impression/ reports	N	Y
213	WGA PGD II		Mol- PGD	amplification	Blastomere/ blastocyst	>4 Blastomeres	2-8 °C	2	medical impression/ reports	N	Y
214	WT1 common somatic mutation	AML, ALL, Prostate cancer, Breast cancer	Mol- Trans	Sanger Sequencing	Tis or DNA or FFPE or PBE	5-10 ml	2-8 °C	14	medical compression/ reports/ pathology results/ doctor's prescription	N	Y
215	Y-Microdeletion	Y-Microdeletion; Male Infertility , Azoospermia	Mol- Others	PCR	PBE or DNA	5-10 ml	2-8 °C	14	medical impression/ reports	Y	Y



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N



قابلیت اورژانس

دارد؟

**Does it
have an
emergenc
function**

N

Y



قابلیت اورژانس
دارد؟
**Does it
have an
emergenc
function**

N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
Y
*Y
N
Y
Y



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
N
N
N
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
N
Y
Y
Y
Y
Y
**Y



قابلیت اورژانس دارد؟ Does it have an emergenc function
Y
**Y
N
Y
Y
Y
Y
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
*Y
N
N
N
*Y
N



قابلیت اورژانس دارد؟ Does it have an emergenc function
Y
Y
N
N